

Epinephrine Treatments

2024

AllergyAsthmaNetwork.org • 800-878-4403

Brand Devices

Generic Devices

	Auvi-Q®	EpiPen®	neffy®	Auto-Injector (Amneal)	Auto-Injector (Teva)	Auto-Injector (Viartis)
Type	Auto-Injector	Auto-Injector	Nasal spray	Auto-Injector	Auto-Injector	Auto-Injector
Dosage	0.10 mg for 16.5 - 33 lbs. 0.15 mg for 33 - 66 lbs. 0.3 mg for over 66 lbs.	0.15 mg for 33 - 66 lbs. 0.3 mg for over 66 lbs.	2 mg for over 66 lbs.	0.15 mg for 33 - 66 lbs. 0.3 mg for over 66 lbs.	0.15 mg for 33 - 66 lbs. 0.3 mg for over 66 lbs.	0.15 mg for 33 - 66 lbs. 0.3 mg for over 66 lbs.
Storage Temperature	68 to 77 degrees F	68 to 77 degrees F	68 to 77 degrees F (with excursions up to 122 degrees F)	68 to 77 degrees F	68 to 77 degrees F	68 to 77 degrees F
Administration	Outer middle of thigh	Outer middle of thigh	Spray in nostril (use same nostril if second dose is needed)	Outer middle of thigh	Outer middle of thigh	Outer middle of thigh
Hold Time	2 seconds	3 seconds	None	10 seconds	3 seconds	3 seconds
Is a trainer available?	Yes	Yes	Yes	Yes	Yes	Yes
Twin-packs available?	Yes	Yes	Yes	Yes	Yes	Yes
Shelf life	12 to 18 months	12 to 18 months	30 months	12 to 18 months	12 to 18 months	12 to 18 months
Special feature	Voice prompt; needle fully covered after injection	Needle fully covered after injection	Needle-free	If needle is sticking out after injection, you received the dose	Needle fully covered after injection	Needle fully covered after injection
Manufacturer	Kaléo	Viartis	ARS Pharmaceuticals	Amneal Pharmaceuticals	Teva Pharmaceuticals	Viartis
Website	auvi-q.com	epipen.com	neffy.com	epinephrineautoinject.com	tevaepinephrine.com	epipen.com
Patient assistance	877-302-8847 	800-395-3376 	877-696-3339 	800-934-6729 	844-248-7949 	800-395-3376 